Performanc

FOXP3 Monoclonal Antibody (FJK-16s), eFluor™ 450, eBioscience™

Size	100 µg
Species Reactivity	Bovine, Dog, Cat, Mouse, Pig, Rat
Published Species	Rat, Bacteria, Mouse, Human
Host/Isotype	Rat / IgG2a, kappa
Recommended Isotype Control	Rat IgG2a kappa Isotype Control (eBR2a), eFluor™ 450, eBioscience™
Class	Monoclonal
Туре	Antibody
Clone	FJK-16s
Conjugate	eFluor™ 450
Excitation/Emission Max	405/445 nm
Form	Liquid
Concentration	0.2 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_1518812

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	0.06 µg/test	96 Publications
ELISA (ELISA)	-	1 Publication

Product Specific Information

Description: The FJK-16s antibody reacts with mouse, rat, dog, porcine, bovine and cat Foxp3 also known as FORKHEAD BOX P3, SCURFIN, and JM2; cross reactivity of this antibody to other proteins has not been determined. Foxp3, a 49-55 kDa protein, is a member of the forkhead/winged-helix family of transcriptional regulators, and was identified as the gene defective in 'scurfy' (sf) mice. Constitutive high expression of foxP3 mRNA has been shown in CD4+CD25+ regulatory T cells (Treg cells), and ectopic expression of foxp3 in CD4+CD25- cells imparts a Treg phenotype in these cells.

Immunoblotting with FJK-16s antibody has mapped the epitope to amino acids 75-125 of the mouse Foxp3 protein. In the human, this region has been shown to be alternatively spliced at the mRNA level. Both the alternatively-spliced and non-spliced isoforms are present in the CD4+CD25+ subset of lymphocytes. Preliminary RT-PCR experiments have not revealed this alternatively-spliced isoform in mouse splenocytes, suggesting different gene regulation in the mouse and human.

Please note that FJK-16s has been optimized for use with the Foxp3 Transcription Factor Buffer Staining Set (Product # 00-5523-00). The use of other fixation and staining buffers is not recommended.

Applications Reported: This FJK-16s antibody has been reported for use in intracellular staining followed by flow cytometric analysis.

Applications Tested: This FJK-16s antibody has been tested by intracellular staining using Foxp3 buffers of mouse splenocytesusing Foxp3/Transcription Factor Buffer Set (Product # 00-5523-00) and protocol. Please see Best Protocols Section (Staining intracellular Antigens for Flow Cytometry) for staining protocol (refer to Protocol B: One-step protocol for intracellular (nuclear) proteins). This can be used at less than or equal to 0.06 μ g per test. A test is defined as the amount (μ g) of antibody that will stain a cell sample in a final volume of 100 μ L. Cell number should be determined empirically but can range from 10^5 to 10^8 cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.

eFluor® 450 is an alternative to Pacific Blue®. eFluor® 450 emits at 445 nm and is excited with the Violet laser (405 nm). Please make sure that your instrument is capable of detecting this fluorochome.

Excitation: 405 nm; Emission: 445 nm; Laser: Violet Laser.

Filtration: 0.2 µm post-manufacturing filtered.

Product Images For FOXP3 Monoclonal Antibody (FJK-16s), eFluor™ 450, eBioscience™



FOXP3 Antibody (48-5773-82)

Intracellular staining of mouse splenocytes. As expected based on known relative expression patterns, Foxp3 clone FJK-16s stains a subset of the CD4+ T cells and does not stain the CD8+ T cells. Details: Balb/c splenocytes were surface stained with CD3 (clone 17A2), CD4 (clone GK1.5) and CD8 (clone 53-6.7), followed by intracellular staining with Foxp3 (clone FJK-16s) using the Foxp3 /Transcription Factor Staining Buffer Set and protocol. Lymphocytes in the CD3+CD8+ (blue histogram) and CD3+CD4+ (purple histogram) gates were used for analysis. {RE}

FOXP3 Antibody (48-5773-82) in Flow

Surface staining of C57BL/6 splenocytes with Anti-Mouse CD4 FITC (Product # 11-0041-82) followed by fixation and permeabilization with the Foxp3 Staining Buffers (Product # 00-5523-00) and intracellular staining with 0.03 µg of Rat IgG2a kappa Isotype Control eFluor® 450 (Product # 48-4321-82) (left) or 0.03 µg of Anti-Mouse/Rat Foxp3 eFluor® 450 (right). Cells in the lymphocyte gate were used for analysis.

97 References

Flow Cytometry (96)

iScience Ischemic Stroke Impacts the Gut Microbiome, Ileal Epithelial and Immune Homeostasis. "48-5773-82 was used in Flow cytometry/Cell sorting to provide a comprehensive overview of stroke-induced intestinal dysbiosis and unveil responsive gene programming within iECs with implications for disease development." Authors: Ge Y,Zadeh M,Yang C,Candelario-Jalil E,Mohamadzadeh M	
Cell reports Engineered human cytokine/antibody fusion proteins expand regulatory T cells and confer autoimmune disease protection. "48-5773-82 was used in Flow cytometry/Cell sorting to provide a roadmap for IC design and establish a Treg-biased immunotherapy that could be clinically translated for autoimmune disease treatment." Authors: VanDyke D,Iglesias M,Tomala J,Young A,Smith J,Perry JA,Gebara E,Cross AR,Cheung LS,Dykema AG, Orcutt-Jahns BT,Henclová T,Golias J,Balolong J,Tomasovic LM,Funda D,Meyer AS,Pardoll DM,Hester J,Issa F,Hunter CA,Anderson MS,Bluestone JA,Raimondi G,Spangler JB	Year 2022 Species Mouse Dilution 1:50

View more Flow references on thermofisher.cn

ELISA (1)

Journal of immunotherapy (Hagerstown, Md. : 1997) High-avidity autoreactive CD4+ T cells induce host CTL, overcome T	Year 2009
(regs) and mediate tumor destruction.	
"48-5773 was used in an ELISpot assay to develop a model of melanoma immunotherapy focused on the exploitation of high-avidity CD4+ T cells that become generated in germline antigen-deficient mice."	Mouse
Authors: Brandmaier AG,Leitner WW,Ha SP,Sidney J,Restifo NP,Touloukian CE	

More applications with references on thermofisher.cn

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. No OTHER WARRANTES, EXPRESS or RIMPLED, ARRANTES, EXPRESS TOR NAVP PARTICULAR PRENESS or RIMPLED, ARRANTES PORTES VERCHANTABILITY. THRESS FOR ANY PARTICULAR PRENESS or RIMPLED, ARRANTES PERFESS OR REFUND FOR THE NON-CONFORMING PRODUCT(s) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCT(s) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OF OR BY BUYER. (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HADDLING OF THE PRODUCTS. UNder Set on the Product or in the documentation accompanying the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to numar or animals.