Performance guarenteed'


AMPD2 Polyclonal Antibody

Product Details

| Size | 400 μL |
|--------------------|---|
| Species Reactivity | Human |
| Published Species | Human |
| Host/Isotype | Rabbit / IgG |
| Class | Polyclonal |
| Туре | Antibody |
| Conjugate | Unconjugated |
| Immunogen | KLH conjugated synthetic peptide between 187-217 amino acids from the central region of human AMPD2 |
| Form | Liquid |
| Concentration | 2.0 mg/mL |
| Purification | Ammonium sulfate precipitation, Size-exclusion - Dialysis |
| Storage buffer | PBS |
| Contains | 0.09% sodium azide |
| Storage conditions | -20° C, Avoid Freeze/Thaw Cycles |
| RRID | AB_2543627 |
| | |

| Applications | Tested Dilution | Publications |
|---|-----------------|---------------|
| Western Blot (WB) | 1:1,000 | 1 Publication |
| Immunohistochemistry (Paraffin) (IHC (P)) | 1:50-1:100 | - |
| Flow Cytometry (Flow) | 1:10-1:50 | 1 Publication |
| Immunoprecipitation (IP) | - | 1 Publication |


Product Images For AMPD2 Polyclonal Antibody



AMPD2 Antibody (PA5-26127) in IHC (P)


Immunohistochemical analysis of formalin-fixed, paraffin-embedded human skeletal muscle tissue using an AMPD2 polyclonal antibody (Product # PA5-26127), followed by HRP-conjugated secondary antibody and DAB staining.

1



AMPD2 Antibody (PA5-26127) in WB

Western blot analysis using an AMPD2 polyclonal antibody (Product # PA5-26127) in K562 cell lysates (35 µg per lane).



AMPD2 Antibody (PA5-26127) in Flow

Flow cytometry analysis of K562 cells using an AMPD2 polyclonal antibody (Product # PA5-26127) (bottom) compared to a negative control cell (top) at a dilution of 1:10-50, followed by a FITC-conjugated goat anti-rabbit antibody

View more figures on thermofisher.cn

3 References

Western Blot (1)

International journal of molecular sciences

The Anti-Glucocorticoid Receptor Antibody Clone 5E4: Raising Awareness of Unspecific Antibody Binding.

"Published figure using AMPD2 polyclonal antibody (Product # PA5-26127) in Western Blot" Authors: Ehlers L,Kirchner M,Mertins P,Strehl C,Buttgereit F,Gaber T

Flow Cytometry (1)

| nmunoprecipitation (1) | |
|--|-------------------------|
| Authors: Ehlers L,Kuppe A,Damerau A,Wilantri S,Kirchner M,Mertins P,Strehl C,Buttgereit F,Gaber T | |
| "PA5-26127 was used in Immunoprecipitation, Flow Cytometry to identify eAMPD2 as a novel regulator of the extracellular ATP-adenosine balance adding to the immunomodulatory CD39-CD73 system." | |
| adenine nucleotide metabolism. | Species Human |
| Surface AMP deaminase 2 as a novel regulator modifying extracellular | Species |
| FASEB journal : official publication of the Federation of American Societies for Experimental Biology | Year 2021 |

| FASEB journal : official publication of the Federation of American | Year |
|--|----------------|
| Societies for Experimental Biology | 2021 |
| Surface AMP deaminase 2 as a novel regulator modifying extracellular | Species |
| adenine nucleotide metabolism. | Human |
| "PA5-26127 was used in Immunoprecipitation, Flow Cytometry to identify eAMPD2 as a novel regulator of the extracellular ATP-adenosine balance adding to the immunomodulatory CD39-CD73 system." Authors: Ehlers L,Kuppe A,Damerau A,Wilantri S,Kirchner M,Mertins P,Strehl C,Buttgereit F,Gaber T | |

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation For research use uniquestic procedures. Not for less authorization, "Policial set warranty expedited by Porton substrated to perfavor operator operator operator operator operator in publications in effect at the time of sauce set ("bottom in the Production and use as for the interview and the production and use as for the interview and the production and use as for the product set warranty is good and does not expense in the any Product will conform to such all the product set warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. No OTHER WARRANTES, EXPRESS OR IMPED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTES OF MERCHANTABILITY, ITHESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS) ON LINE TO FORCE MAJELER (") MISUE, FAULT OR NEGLIGENCE OF OR BY DURCT, (ii) USE OF THE PRODUCTS IN A MANEER FOR WHICH THEY WERE NOT DESIGNED, OR ((ii) MIROPER STORAGE AND HARD/CONFORMING PRODUCTS). Unless otherwise expressly stated on the Product or in the documentation accompanying the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro detarrent of a companying the product is intended for research only and is not to be used for any other purpose, including without limitation, and and accompanying the product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro detarrent of a companying and accompanying the product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro detarrent of a vide as a limited to hus diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals

Year 2022