Performanc guarenteed

# CD64 Monoclonal Antibody (X54-5/7.1), APC-eFluor™ 780, eBioscience™

Product Details	
Size	100 µg
Species Reactivity	Mouse
Host/Isotype	Mouse / IgG1, kappa
Recommended Isotype Control	Mouse IgG1 kappa Isotype Control (P3.6.2.8.1), APC-eFluor™ 780, eBioscience™
Class	Monoclonal
Туре	Antibody
Clone	X54-5/7.1
Conjugate	APC-eFluor™ 780
Excitation/Emission Max	756/785 nm
Immunogen	BALB/c mouse CD64-human IgG Fc fusion protein
Form	Liquid
Concentration	0.2 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_2735012

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	0.5 µg/test	1 Publication

## **Product Specific Information**

Description: This X54-5/7.1 monoclonal antibody recognizes allelic variants a and b of mouse CD64, a high affinity Fc binding receptor also known as Fc-gammaR1 or Fcgr1. Alloform a is present in the following mouse strains: C57BL/6, BALB/c, DBA/2, and alloform b in strains: C3H/HeJ, CBA/J, NZW, SJL/J, 129/SvJ. In addition this antibody has been reported positive in strains: AKR, ALR, BUB, C58, CE, HRS,MRL, MON, NZB, NZO, NZW, PL, SJL, ST and SWR. This clone X54-5/7.1 will not recognize the allelic variant d, present in NOD/Lt mice. This antibody has been reported to stain RAW264.7 cell line.

Applications Reported: This X54-5/7.1 antibody has been reported for use in flow cytometric analysis.

Applications Tested: This X54-5/7.1 antibody has been tested by flow cytometric analysis of mouse bone marrow cells. This may be used at less than or equal to 0.5  $\mu$ g per test. A test is defined as the amount ( $\mu$ g) of antibody that will stain a cell sample in a final volume of 100  $\mu$ L. Cell number should be determined empirically but can range from 10^5 to 10^8 cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.

APC-eFluor<sup>™</sup> 780 emits at 780 nm and is excited with the Red laser (633 nm). Please make sure that your instrument is capable of detecting this fluorochrome.


Light sensitivity: This tandem dye is sensitive to photo-induced oxidation. Please protect this vial and stained samples from light.

1

Fixation: Samples can be stored in IC Fixation Buffer (Product # 00-8222) (100 µL of cell sample + 100 µL of IC Fixation Buffer) or 1-step Fix/Lyse Solution (Product # 00-5333) for up to 3 days in the dark at 4°C with minimal impact on brightness and FRET efficiency/compensation. Some generalizations regarding fluorophore performance after fixation can be made, but clone specific performance should be determined empirically.


Excitation: 633-647 nm; Emission: 780 nm; Laser: Red Laser

# Product Images For CD64 Monoclonal Antibody (X54-5/7.1), APC-eFluor™ 780, eBioscience™



#### CD64 Antibody (47-0641-82)

Staining of mouse bone marrow. As expected based on the known relative expression patterns, CD64 clone X54-5/7.1 stains monocytes with the highest intensity (purple histogram) and B cell precursors with lower intensity (orange and blue histograms). Details: Mouse bone marrow was Fc blocked, stained with CD64 (clone X54-5/7.1), and co-stained with CD45R (B220) (clone RA3-6B2) and CD43 (clone eBioR2/60). Cells in the indicated color gates were used for analysis. {RE}



### CD64 Antibody (47-0641-82) in Flow

BALB/c mouse bone marrow cells were Fc blocked with CD16/CD32 Antibody (Product # 14-0161-82) and normal mouse serum, and stained with CD11b Monoclonal Antibody, FITC (Product # 11-0112-82) and 0.25 µg of Mouse IgG1 kappa Isotype Control, APC-eFluor 780 (Product # 47-4714-82) (left) or 0.25 µg of CD64 Monoclonal Antibody, APC-eFluor 780 (right). All cells were used for analysis.

# □ 1 Reference

## Flow Cytometry (1)

Science advances	Year
Circulating hemopexin modulates anthracycline cardiac toxicity in patients and in mice.	2022
"Published figure using CD64 monoclonal antibody (Product # 47-0641-82) in Flow Cytometry"	
Authors: Liu J,Lane S,Lall R,Russo M,Farrell L,Debreli Coskun M,Curtin C,Araujo-Gutierrez R,Scherrer-Crosbie M, Trachtenberg BH.Kim J.Tolosano E.Ghigo A.Gerszten RE.Asnani A	

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of said, as set forth in the Production documentation by a specifications and effect at the time of said said. The set of the set of substantially is a set of the product specifications in effect at the time of said. The product specifications in effect at the time of said. The product specifications is made. The warranty provide therein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty provide therein is valid only only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty provide therein is valid only only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty provide therein is valid only only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty provide therein is valid only only when used by properly trained individuals. Unless otherwise stated on anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. No OTHER WARRANTES, EXPRESS OR INNELD, ARE ORANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTES OF ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEMPT POR NON-CONFORMING PRODUCTS DURING THE WARRANTE PROVIDENT EXPRESS OR INFRINGEMENT. BUYER'S EXCLUSIVE REMEMPT POR NOVE CONFORMING PRODUCTS ON STATE, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS) AT SETLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (III) MISUE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A

HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the docum ompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized com diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals