




SELENBP1 Monoclonal Antibody (OTI3H1)

| Product Details | |
|------------------------|---------------------------------------------------------------------------|
| Size | 100 μL |
| Species Reactivity | Human, Mouse |
| Host/Isotype | Mouse / IgG2a |
| Class | Monoclonal |
| Туре | Antibody |
| Clone | OTI3H1 |
| Conjugate | Unconjugated |
| Immunogen | Full length human recombinant protein of SELENB1 produced in HEK293T cell |
| Form | Liquid |
| Concentration | 1 mg/mL |
| Purification | Affinity Chromatography |
| Storage buffer | PBS, pH 7.3, with 1% BSA, 50% glycerol |
| Contains | 0.02% sodium azide |
| Storage conditions | -20° C, Avoid Freeze/Thaw Cycles |
| RRID | AB_2723626 |


| Applications | Tested Dilution | Publications |
|-------------------------------------------|-----------------|--------------|
| Western Blot (WB) | 1:2,000 | - |
| Immunohistochemistry (Paraffin) (IHC (P)) | 1:150 | - |

Product Images For SELENBP1 Monoclonal Antibody (OTI3H1)




SELENBP1 Antibody (MA5-25953) in IHC (P)

Immunohistochemistry was performed on paraffin-embedded adenocarcinoma of human colon tissue. To expose target proteins, 10mM citric buffer, pH6.0, 120°C for 3min was used. Following antigen retrieval, tissues were probed with a SELENBP1 monoclonal antibody (Product # MA5-25953).



SELENBP1 Antibody (MA5-25953) in IHC (P)

Immunohistochemistry was performed on paraffin-embedded carcinoma of human kidney tissue. To expose target proteins, 10mM citric buffer, pH6.0, 120° C for 3min was used. Following antigen retrieval, tissues were probed with a SELENBP1 monoclonal antibody (Product # MA5-25953).



SELENBP1 Antibody (MA5-25953) in IHC (P)

Immunohistochemistry was performed on paraffin-embedded human lymph node tissue. To expose target proteins, 10mM citric buffer, pH6.0, 120°C for 3min was used. Following antigen retrieval, tissues were probed with a SELENBP1 monoclonal antibody (Product # MA5-25953).

View more figures on thermofisher.cn

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization, Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless of service in the product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICUAR PURPOSE, OR NON INFRINGEMENT.

BUYER'S EXCLUSIVE REMEDY FOR NON-CORPORIMO PRODUCTS DURING THE WARRANTY PERFIOI DIS LIMITED, ARE CREATED INCLUDING WITHOUT LIMITATION, IMPLED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICUARD PURPOSE, OR NON INFRINGEMENT.

BUYER'S EXCLUSIVE REMEDY FOR NON-CORPORIMO PRODUCTS DURING THE WARRANTY PERFIOI DIS LIMITED, ARE SUBJECT OF REFUND FOR THE NON-CONFORMING PRODUCTS, OR AT SUBJECT SOLE OPTION. THERE IS NO ADMINISTRATION TO REPLAIR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS, OR AT SUBJECT SOLE OPTION. THERE IS NO ADMINISTRATION TO SUBJECT SOLE OPTION.

HANDLING OF THE PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any ot