

CD274 (PD-L1, B7-H1) Monoclonal Antibody (MIH1), PE-Cyanine5, eBioscience™

Product Details	
Size	100 Tests
Species Reactivity	Human
Host/Isotype	Mouse / IgG1, kappa
Recommended Isotype Control	Mouse IgG1 kappa Isotype Control (P3.6.2.8.1), PE-Cyanine5, eBioscience™
Class	Monoclonal
Type	Antibody
Clone	MIH1
Conjugate	PE-Cyanine5
Excitation/Emission Max	568/666 nm
Form	Liquid
Concentration	5 µL/Test
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2, with 0.2% BSA
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_2802211

Applications	Tested Dilution	Publications
Western Blot (WB)	-	1 Publication
Immunohistochemistry (IHC)	-	1 Publication
Immunohistochemistry (Paraffin) (IHC (P))	-	1 Publication
Flow Cytometry (Flow)	5 µL (0.25 µg)/test	23 Publications

Product Specific Information

Description: The MIH1 monoclonal antibody reacts with human B7-H1, also known as PD-L1. B7-H1, a member of the B7 family, has a predicted molecular weight of approximately 40 kDa and belongs to the Ig superfamily. B7-H1 is expressed on a majority of leukocytes. B7-H1 is a ligand for PD-1. Interaction of PD-1 with either PD-L1 (B7-H1) or PD-L2 (B7-DC) results in inhibition of T and B cell responses. MIH1 is reported to be a blocking antibody.

Applications Reported: This MIH1 antibody has been reported for use in flow cytometric analysis.

Applications Tested: This MIH1 antibody has been pre-diluted and tested by flow cytometric analysis of stimulated normal human peripheral blood cells. This may be used at 5 µL (0.25 µg) per test. A test is defined as the amount (µg) of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 10⁵ to 10⁸ cells/test.


Light sensitivity: This tandem dye is sensitive to photo-induced oxidation. Please protect this vial and stained samples from light.

Fixation: Samples can be stored in IC Fixation Buffer (Product # 00-8222-49) (100 µL of cell sample + 100 µL of IC Fixation

Buffer) or 1-step Fix/Lyse Solution (Product # 00-5333-57) for up to 3 days in the dark at 4°C with minimal impact on brightness and FRET efficiency/compensation. Some generalizations regarding fluorophore performance after fixation can be made, but clone specific performance should be determined empirically.

Excitation: 488-561 nm; Emission: 667 nm; Laser: Blue Laser, Green Laser, Yellow-Green Laser

Product Images For CD274 (PD-L1, B7-H1) Monoclonal Antibody (MIH1), PE-Cyanine5, eBioscience™



CD274 (PD-L1, B7-H1) Antibody (15-5983-42) in Flow
Normal human peripheral blood cells were unstimulated (left) or stimulated with PHA-L overnight (right). Cells were then stained with Mouse IgG1 kappa Isotype Control, PE-Cyanine5 (Product # 15-4714-81) (blue histogram) or CD274 Monoclonal Antibody, PE-Cyanine5 (purple histogram). Viable cells in the lymphocyte gate were used for analysis, as determined by Fixable Viability Dye eFluor 780 (Product # 65-0865-18).

[View more figures on thermofisher.cn](https://thermofisher.cn)

Western Blot (1)

<p>Cancers</p> <p>Microbe-Mediated Activation of Toll-like Receptor 2 Drives PDL1 Expression in HNSCC.</p> <p>"Published figure using CD274 (PD-L1, B7-H1) monoclonal antibody (Product # 15-5983-42) in Western Blot"</p> <p>Authors: Mann JE,Ludwig ML,Kulkarni A,Scheftz EB,Murray IR,Zhai J,Gensterblum-Miller E,Jiang H,Brenner JC</p>	<p>Year 2021</p>
---	-----------------------------

Immunohistochemistry (1)

<p>Biomedicines</p> <p>Concomitant High Apoptosis Inhibitor of Macrophage (AIM) and Low Prostate-Specific Antigen (PSA) Indicates Activated T Cell-Mediated Anticancer Immunity, Enhance Sensitivity to Pembrolizumab, and Elicit Good Prognosis in Prostate Cancer.</p> <p>"Published figure using CD274 (PD-L1, B7-H1) monoclonal antibody (Product # 15-5983-42) in Immunohistochemistry"</p> <p>Authors: Bamodu OA,Wang YH,Yeh CT,Ho CH,Chiang YT,Kao WT,Liu CH,Wu CC</p>	<p>Year 2021</p>
---	-----------------------------

Immunohistochemistry (Paraffin) (1)

<p>International journal of molecular sciences</p> <p>CD24: A Marker for an Extended Expansion Potential of Urothelial Cancer Cell Organoids In Vitro?</p> <p>"Published figure using CD274 (PD-L1, B7-H1) monoclonal antibody (Product # 15-5983-42) in Immunohistochemistry (Paraffin)"</p> <p>Authors: Geng R,Harland N,Montes-Mojarro IA,Fend F,Aicher WK,Stenzl A,Amend B</p>	<p>Year 2022</p>
--	-----------------------------

Flow Cytometry (23)

<p>Frontiers in immunology</p> <p>GM-CSF, Flt3-L and IL-4 affect viability and function of conventional dendritic cell types 1 and 2.</p> <p>"Published figure using CD274 (PD-L1, B7-H1) monoclonal antibody (Product # 15-5983-42) in Flow Cytometry"</p> <p>Authors: Lellahi SM,Azeem W,Hua Y,Gabriel B,Paulsen Rye K,Reikvam H,Kalland KH</p>	<p>Year 2023</p>
---	-----------------------------

[View more Flow references on thermofisher.cn](#)

More applications with references on thermofisher.cn

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (i) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (ii) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (iii) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (iv) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.