

IL-8 (CXCL8) Monoclonal Antibody (790A28G2)

Product Details	
Size	500 µg
Species Reactivity	Human
Published Species	Human
Host/Isotype	Mouse / IgG1, kappa
Class	Monoclonal
Type	Antibody
Clone	790A28G2
Conjugate	Unconjugated
Immunogen	Purified recombinant human IL-8.
Form	Liquid
Concentration	1 mg/mL
Purification	Protein A
Storage buffer	PBS, pH 7.4, with 5mg/mL BSA
Contains	0.5mg/mL sodium azide
Storage conditions	4° C
RRID	AB_2536260

Applications	Tested Dilution	Publications
ELISA (ELISA)	Assay-dependent	6 Publications

Product Specific Information

This product is recommended as a detection antibody in Sandwich ELISA applications.

This clone can also be referred to as clone number 790A 28G2 1E3.

ELISA (6)

<p>International journal of molecular sciences</p> <p>IL-8 and cathepsin B as melanoma serum biomarkers.</p> <p>"AHC0784 was used in ELISA to study IL8 and Cathepsin B in the prognosis of melanoma mortality."</p> <p>Authors: Zhang H,Fu T,McGettigan S,Kumar S,Liu S,Speicher D,Schuchter L,Xu X</p>	<p>Year 2015</p> <p>Species Human</p>
<p>Journal of leukocyte biology</p> <p>The role of the Src family kinase Lyn in the immunomodulatory activities of cathelicidin peptide LL-37 on monocytic cells.</p> <p>"AHC0784 was used in ELISA to report that SFKs regulate the chemokine induction activity of LL-37 in monocytic cells."</p> <p>Authors: Nijnik A,Pistolici J,Cho P,Filewod NC,Falsafi R,Ramin A,Harder KW,Hancock RE</p>	<p>Year 2012</p> <p>Species Human</p>

[View more ELISA references on thermofisher.cn](#)

More applications with references on thermofisher.cn

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.