

FAS Monoclonal Antibody (DX2)

Product Details

Size	50 µg
Species Reactivity	Human
Published Species	Human
Host/Isotype	Mouse / IgG1
Class	Monoclonal
Type	Antibody
Clone	DX2
Conjugate	Unconjugated
Immunogen	Human Fas/CD95 cDNA.
Form	Liquid
Concentration	1.0 mg/mL
Purification	Protein A
Storage buffer	0.01M PBS, pH 7.4
Contains	15mM sodium azide
Storage conditions	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.
RRID	AB_558307

Applications

Tested Dilution

Publications

Flow Cytometry (Flow)	4-20 µg/mL	6 Publications
Miscellaneous PubMed (Misc)	-	4 Publications

Product Specific Information

Recommended positive controls: Raji.

Store product as a concentrated solution. Centrifuge briefly prior to opening the vial.

Flow Cytometry (6)

The American journal of pathology

Characterization and Purification of Neoplastic Cells of Nodular Lymphocyte Predominant Hodgkin Lymphoma from Lymph Nodes by Flow Cytometry and Flow Cytometric Cell Sorting.

"MA120163 was used in flow cytometry to document and describe lymphocyte predominant cells from lymph nodes involved in nodular lymphocyte predominant Hodgkin lymphoma"

Authors: Fromm JR,Thomas A,Wood BL

Year
2017

Species
Human

Blood

Flow sorting and exome sequencing reveal the oncogenome of primary Hodgkin and Reed-Sternberg cells.

"MA1-20163 was used in flow cytometry to identify alterations in genes involved in antigen presentation, chromosome integrity, transcriptional regulation, and ubiquitination in Hodgkin and Reed-Sternberg cells"

Authors: Reichel J,Chadburn A,Rubinstein PG,Giulino-Roth L,Tam W,Liu Y,Gaiolla R,Eng K,Brody J,Inghirami G,Carlo-Stella C,Santoro A,Rahal D,Totonchy J,Elemento O,Cesarman E,Roshal M

Year
2015

Species
Human

[View more Flow references on thermofisher.cn](#)

Miscellaneous PubMed (4)

Human immunology

Increased caspase activity primes human Lyme arthritis synovial T cells for proliferation and death.

"MA1-20163 was used in flow cytometry to study the turnover dynamics of human T cells."

Authors: Thai PT,Collins CC,Fortner KA,Koenig A,Hayes SM,Budd RC

Year
2011

Species
Human

Journal of immunology (Baltimore, Md. : 1950)

Early induction of polyfunctional simian immunodeficiency virus (SIV)-specific T lymphocytes and rapid disappearance of SIV from lymph nodes of sooty mangabeys during primary infection.

"MA1-20163 was used in flow cytometry to study the contribution of the cellular immune response in maintaining nonpathogenic SIV infection in sooty mangabeys."

Authors: Meythaler M,Wang Z,Martinot A,Pryputniewicz S,Kasheta M,McClure HM,O'Neil SP,Kaur A

Year
2011

Species
Human

[View more Misc references on thermofisher.cn](#)

More applications with references on thermofisher.cn

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.