


SNRPB Monoclonal Antibody (Y12)

Product Details

Size	500 µL
Species Reactivity	Fruit fly, Human, Mouse, Pig, Rabbit, Rat
Published Species	Dog, Rat, Fruit fly, Hamster, Human, Mouse, Flatworm
Host/Isotype	Mouse / IgG3, kappa
Class	Monoclonal
Type	Antibody
Clone	Y12
Conjugate	Unconjugated
Form	Liquid
Concentration	200 mg/mL
Purification	Protein G
Storage buffer	PBS, pH 7.4, with 0.2% BSA
Contains	0.09% sodium azide
Storage conditions	4° C
RRID	AB_10944191


Applications	Tested Dilution	Publications
Western Blot (WB)	0.5-1.0 µg/mL	5 Publications
Immunohistochemistry (IHC)	-	4 Publications
Immunohistochemistry (Paraffin) (IHC (P))	1:50-1:100	-
Immunohistochemistry (Frozen) (IHC (F))	-	1 Publication
Immunocytochemistry (ICC/IF)	2 µg/mL	7 Publications
Immunoprecipitation (IP)	-	7 Publications
in situ PLA (PLA)	-	1 Publication
Immunomicroscopy (IM)	Assay-dependent	-
RNA Immunoprecipitation (RIP)	1.5 µg/10 ⁶ cells	-

Product Images For SNRPB Monoclonal Antibody (Y12)




SNRPB Antibody (MA5-13449)

Antibody specificity was demonstrated by siRNA mediated knockdown of target protein. HeLa cells were transfected with SNRPB siRNA and reduction of specific signal was observed in immunofluorescence application using SNRPB Mouse Monoclonal antibody (Product # MA5-13449). {KD}



SNRPB Antibody (MA5-13449) in WB

Western blot analysis was performed on Nuclear enriched extracts (30 µg lysate) of HeLa (Lane 1), NTERA-2 (Lane 2), NIH/3T3 (Lane 3), Hep G2 (Lane 4), K-562 (Lane 5), SH-SY5Y (Lane 6), HEK-293 (Lane 7) and A-431 (Lane 8). The blots were probed with Mouse Anti-SNRPB Monoclonal Antibody (Product # MA5-13449, 1 µg/mL) and detected by chemiluminescence using Goat anti-Mouse IgG (H+L) Superclonal™ Secondary Antibody, HRP conjugate (Product # A28177, 0.25 µg/mL, 1:4000 dilution). A 24 kDa and 30 kDa band corresponding to SNRPB was observed across the cell lines tested. Known quantity of protein samples were electrophoresed using Novex® NuPAGE® 4-12 % Bis-Tris gel (Product # NP0321BOX), XCell SureLock™ Electrophoresis System (Product # EI0002) and Novex® Sharp Pre-Stained Protein Standard (Product # LC5800). Resolved proteins were then transferred onto a nitrocellulose membrane with iBlot® 2 Dry Blotting System (Product # IB21001). The membrane was probed with the relevant primary antibody following blocking with 5 % skimmed milk. Chemiluminescent detection was performed using Pierce™ ECL Western Blotting Substrate (Product # 32106).



SNRPB Antibody (MA5-13449)

Antibody specificity was demonstrated by siRNA mediated knockdown of target protein. HeLa cells were transfected with SNRPB siRNA and decrease in signal intensity was observed in western blot application using anti-SNRPB siRNA (Product # MA5-13449). {KD}

View more figures on thermofisher.cn

Western Blot (5)

<p>Nature communications</p> <p>Therapeutic manipulation of IKBKAP mis-splicing with a small molecule to cure familial dysautonomia.</p> <p>"MA5-13449 was used in Immunocytochemistry to indicate that small synthetic molecules affecting phosphorylation state of SRSFs is available as a new therapeutic modality for mechanism-oriented precision medicine of splicing diseases."</p> <p>Authors: Ajiro M, Awaya T, Kim YJ, Iida K, Denawa M, Tanaka N, Kurosawa R, Matsushima S, Shibata S, Sakamoto T, Studer L, Krainer AR, Hagiwara M</p>	<p>Year 2021</p> <p>Species Human</p> <p>Dilution 1:500</p>
<p>Cell reports</p> <p>Tau-Mediated Disruption of the Spliceosome Triggers Cryptic RNA Splicing and Neurodegeneration in Alzheimer's Disease.</p> <p>"Published figure using SNRPB monoclonal antibody (Product # MA5-13449) in Immunohistochemistry"</p> <p>Authors: Hsieh YC, Guo C, Yalamanchili HK, Abreha M, Al-Ouran R, Li Y, Dammer EB, Lah JJ, Levey AI, Bennett DA, De Jager PL, Seyfried NT, Liu Z, Shulman JM</p>	<p>Year 2019</p>

[View more WB references on thermofisher.cn](#)

Immunohistochemistry (4)

<p>Cell reports</p> <p>Tau-Mediated Disruption of the Spliceosome Triggers Cryptic RNA Splicing and Neurodegeneration in Alzheimer's Disease.</p> <p>"Published figure using SNRPB monoclonal antibody (Product # MA5-13449) in Immunohistochemistry"</p> <p>Authors: Hsieh YC, Guo C, Yalamanchili HK, Abreha M, Al-Ouran R, Li Y, Dammer EB, Lah JJ, Levey AI, Bennett DA, De Jager PL, Seyfried NT, Liu Z, Shulman JM</p>	<p>Year 2019</p>
<p>Human molecular genetics</p> <p>Mutations in splicing factor PRPF3, causing retinal degeneration, form detrimental aggregates in photoreceptor cells.</p> <p>"MA5-13449 was used in immunohistochemistry to investigate the role of mutations in splicing factor PRPF3 in retinal degeneration of photoreceptor cells"</p> <p>Authors: Comitato A, Spanpanato C, Chakarova C, Sanges D, Bhattacharya SS, Marigo V</p>	<p>Year 2007</p> <p>Species Human</p> <p>Dilution 1:100</p>

[View more IHC references on thermofisher.cn](#)

More applications with references on thermofisher.cn

IHC (F) (1) ICC/IF (7) IP (7) PLA (1)

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.