

SARS-CoV-2 Spike Protein S1 Chimeric Recombinant Human Monoclonal Antibody (H6)


Product Details	
Size	100 µL
Species Reactivity	Virus
Host/Isotype	Human / IgG1, kappa
Expression system	HEK293 cells
Class	Recombinant Monoclonal
Type	Antibody
Clone	H6
Conjugate	Unconjugated
Immunogen	Recombinant Human Novel Coronavirus Spike glycoprotein (S) (16-685aa)
Form	Liquid
Concentration	1 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.4, with 50% glycerol
Contains	0.03% ProClin 300
Storage conditions	-20°C or -80°C if preferred
RRID	AB_2866552

Applications	Tested Dilution	Publications
ELISA (ELISA)	1:10,000-1:50,000	-
Lateral Flow (LF)	1:500-1:25,000	-

Product Specific Information


Recombinant anti-SARS-CoV-2 spike is expressed from 293 cells (HEK293). It is a chimeric monoclonal antibody combining the mouse variable regions with the constant domains of human IgG1.

Product Images For SARS-CoV-2 Spike Protein S1 Chimeric Recombinant Human Monoclonal Antibody (H6)



SARS-CoV-2 Spike Protein S1 Chimeric Antibody (MA5-35940) in ELISA

A functional ELISA was performed to by immobilizing SARS-CoV-2-S1-RBD at 2 µg/mL and using SARS-CoV-2 Spike Protein S1 Chimeric Recombinant Monoclonal Antibody (Product # MA5-35940). The EC50 of the antibody is 13.96 - 16.62 ng/mL.



SARS-CoV-2 Spike Protein S1 Chimeric Antibody (MA5-35940) in ELISA

A functional ELISA was performed to by immobilizing SARS-CoV-2-S1-RBD at 2 µg/mL and using SARS-CoV-2 Spike Protein S1 Chimeric Recombinant Monoclonal Antibody (Product # MA5-35940). The EC50 of the antibody is 13.48- 19.50 ng/mL.



SARS-CoV-2 Spike Protein S1 Chimeric Antibody (MA5-35940) in LF

In the Lateral Flow detection system the background of SARS-CoV-2 Spike Protein S1 Chimeric monoclonal antibody (Product # MA5-35940) is clean. The detection limit can be as low as 13.94 ng/mL (0.976 ng/0.07mL), and the sensitivity is exceptional.

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.