



Donkey anti-Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor™ 594

| Product Details | |
|----------------------------|--|
| Size | 1 mg |
| Species Reactivity | Mouse |
| Host/Isotype | Donkey / IgG |
| Class | Polyclonal |
| Туре | Secondary Antibody |
| Conjugate | Alexa Fluor™ 594 |
| Excitation/Emission Max | 590/618 nm |
| Immunogen | Gamma Immunoglobins Heavy and Light chains |
| Form | Liquid |
| Concentration | 2 mg/mL |
| Purification | purified |
| Storage buffer | PBS, pH 7.5 |
| Contains | 5mM sodium azide |
| Storage conditions | 4° C, store in dark |
| RRID | AB_141633 |

| Applications | Tested Dilution | Publications |
|---|-----------------|---------------|
| Western Blot (WB) | - | 0 Publication |
| Immunohistochemistry (IHC) | 1-10 μg/mL | 0 Publication |
| Immunohistochemistry (Paraffin) (IHC (P)) | - | 0 Publication |
| Immunohistochemistry (PFA fixed) (IHC (PFA)) | - | 0 Publication |
| Immunohistochemistry (Frozen) (IHC (F)) | - | 0 Publication |
| Immunohistochemistry - Free Floating (IHC (Free)) | - | 0 Publication |
| Immunocytochemistry (ICC/IF) | 4 μg/mL | 0 Publication |
| Miscellaneous PubMed (Misc) | - | 0 Publication |

Product Specific Information

These donkey anti-mouse IgG (H+L) whole secondary antibodies have been affinity-purified and show minimum crossreactivity to bovine, chicken, goat, quinea pig, hamster, horse, human, rabbit, rat, and sheep serum proteins. Cross-adsorption or pre-adsorption is a purification step to increase specificity of the antibody resulting in higher sensitivity and less background staining. The secondary antibody solution is passed through a column matrix containing immobilized serum proteins from potentially cross-reactive species. Only the nonspecific-binding secondary antibodies are captured in the column, and the highly specific secondaries flow through. The benefits of this extra step are apparent in multiplexing/multicolor-staining experiments (e.g., flow cytometry) where there is potential cross-reactivity with other primary antibodies or in tissue/cell fluorescent staining experiments where there may be the presence of endogenous immunoglobulins.

Alexa Fluor dyes are among the most trusted fluorescent dyes available today. Invitrogen™ Alexa Fluor 594 dye is a bright, red-fluorescent dye with excitation ideally suited to the 594 nm laser line. For stable signal generation in imaging and flow

cytometry, Alexa Fluor 594 dye is pH-insensitive over a wide molar range. Probes with high fluorescence quantum yield and high photostability allow detection of low-abundance biological structures with great sensitivity. Alexa Fluor 594 dye molecules can be attached to proteins at high molar ratios without significant self-quenching, enabling brighter conjugates and more sensitive detection. The degree of labeling for each conjugate is typically 2-8 fluorophore molecules per IgG molecule; the exact degree of labeling is indicated on the certificate of analysis for each product lot.


Using conjugate solutions: Centrifuge the protein conjugate solution briefly in a microcentrifuge before use; add only the supernatant to the experiment. This step will help eliminate any protein aggregates that may have formed during storage, thereby reducing nonspecific background staining. Because staining protocols vary with application, the appropriate dilution of antibody should be determined empirically. For the fluorophore-labeled antibodies a final concentration of 1-10 μ g/mL should be satisfactory for most immunohistochemistry and flow cytometry applications.

Product will be shipped at Room Temperature.

Product Images For Donkey anti-Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor™ 594

Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody (A-21203) in ICC/IF

FGF-2 impedes the AAD-induced anti-EC effect via FGFR1-ERK-MYC signaling. A Cell growth of human ECs receiving the conditioned medium of scrambled- or FGF2 shRNA-transfected NPC tumor cells (n = 5 samples per group). B Representative micrographs of PCNA+ proliferative cells and DAPI signals in ECs treated with vehicle or recombinant human FGF-2. Scale bar = 50 µm. Quantification of PCNA+ signals in mouse and human ECs (n = 8 random fields per group). C Vehicle- or VEGF-treated ECs were challenged with or without sunitinib or FGF-2. Phosphorylation of AKT and ERK in ECs was detected. -actin marks the loading level in each lane (n = 3 samples per group). D QPCR quantification of Fgfr1, Fgfr2, Fgfr3, and Fgfr4 mRNA levels in ECs (n = 3 samples per group). E Vehicle- or FGF-2-treated ECs were challenged with or without various FGFR inhibitors. Phosphorylation of ERK in ECs was detected. actin marks the loading level in each lane (n = 3 samples per group). F Downstream of VEGF signaling transcription factors were selected and detected in vehicle- or FGF-2-treated ECs. Heatmap of qPCR array screened out Myc as the highest upregulated transcription factor. G Correlation of FGF2 and MYC transcriptomic expression levels of human NPCs (NPC, n = 113 samples). Data was extracted from dataset GSE102349. H QPCR quantification of Myc mRNA levels in isolated mouse CD31+ ECs from scrambled- or FGF2 shRNAtransfected NPC tu... Image collected and cropped by CiteAb from the following publication (https://pubmed.ncbi.nlm.nih.gov/35985991), licensed under a CC BY license.




Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody (A-21203) in ICC/IF

Immunofluorescent analysis of MAP2 in the differentiated neurons from H9 ESC-derived NSCs. 2 weeks after differentiation, cells were fixed, permeabilized and stained with a MAP2 rabbit polyclonal antibody (Product # PA5-17646) at 1:100 dilution (green) and a HuC/HuD mouse monoclonal antibody (Product # A-21271), at a concentration of 5 µg/mL (red) in blocking buffer for at least 1 hour at room temperature, and then incubated with goat anti-rabbit IgG secondary antibody, Alexa Fluor Plus 488 conjugate (Product # A32731, green) and a donkey anti-mouse IgG secondary antibody, Alexa Fluor 594 conjugate (Product # A-21203, red) at a dilution of 1:1000 for 1 hour at room temperature. Nuclei (blue) were stained with Hoechst 33342 dye (Product # 62249).

Mouse IgG (H+L) Highly Cross-Adsorbed Secondary Antibody (A-21203) in ICC/IF

Abnormal axonal outgrowth in DM1 hiNeurons at 15 DPI.(a) Representative fluorescent images display axons of control (top panel) and DM1 (bottom panel) hiNeurons visualized by co-staining with TUJ1 (green) and SMI-312 (red) antibodies. Nuclei stained by DAPI. On the right merged images of TUJ1, SMI-312 and DAPI. Scale bar, 50 µm. (b-e) box and whisker plots show measurements of axonal length of control and DM1 hiNeurons. Results were categorized into four categories: the longest/photo (b), medium length/photo (c), the shortest/photo (d) and average/photo (e). Each sample is shown in different color. Line & (+) sign inside the box represent median and mean of replicates (outcome analyzed), respectively. Whiskers show minimum & maximum values. Note that the average length of axons measured in DM1 hiNeurons group is about half length of axons in control group. n = 3 for each group, a total of 36 (the longest, medium and the shortest/photo analysis) or 108 (average/photo analysis) axons were analyzed per sample. A total of 36 images were analyzed per sample. **P<0.01 compared to control group by unpaired t-test. Average axonal length is the average of the longest, medium and the shortest axons per photo. Medium axon is a single representative of the general population of axons excluding the longest and shortest axons in each photo. Image collected and cropped by CiteAb from the following publication (https://pubmed.ncbi.nlm.nih.gov /35776705), licensed under a CC BY license.



View more figures on thermofisher.cn

□ 1545 References

Impact of Micro- and Nano-Plastics on Human Intestinal Organoid-Derived Epithelium. Curr Protoc (2024)

Histological analysis of (antral) follicle density in ovarian cortex tissue attached to stripped endometriomas. J Assist Reprod Genet (2024)

Deficiency of Pdcd10 causes urothelium hypertrophy and vesicle trafficking defects in ureter. FEBS J (2024)

FXTAS Neuropathology Includes Widespread Reactive Astrogliosis and White Matter Specific Astrocyte Degeneration. Ann Neurol (2024)

Melatonin MT1 receptors regulate the Sirt1/Nrf2/Ho-1/Gpx4 pathway to prevent -synuclein-induced ferroptosis in Parkinson's disease. J Pineal Res (2024)

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production obcumentation specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any mode of sample is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample incomplete the product is incompleted to prove the product of the Product or in the documentation accompanying the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, or why or in vito therapeutic uses, or any type of consumption to human or animal diagnostic uses, or why or in vito therapeutic uses, or any type of consumption to human or animal to the buse of the product of the product of in the model and product of the product of the product of in the documentation accompanying the product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, or why or in vito therapeutic uses, or any type of consumption to human or animal to have an animal or an advanced to the product of the prod